



What does it include?

The FitMi System Includes:

- 2 FitMi Pucks
- FitMi Receiver
- Charging Dock
- The RehabStudio App

FlintRehab
Tools to spark recovery

Ready to Start Improving Movement?

To order a FitMi or get more
information, please visit us at:

www.flintrehab.com

You can also contact a
Flint representative at:

Phone: (800) 593-5468

Email: support@flintrehab.com



© 2016 Flint Rehabilitation Devices LLC.

All rights reserved. Reproduction, duplication or modification
in whole or in part is prohibited without prior written permission
from Flint Rehabilitation Devices. All trademarks are owned by or
licensed to Flint Rehabilitation Devices, LLC unless otherwise noted.

Revision 1.28.2016

FitMi
Motion Interface




FlintRehab
Tools to spark recovery



What Apps are Included?

FitMi comes with our unique “RehabStudio” app that turns any computer into a virtual gym.



In RehabStudio, you can create custom therapy regimens from a library of 40 classic exercises. Then, your patients use FitMi to complete the number of repetitions you prescribe for each exercise.

What is FitMi?

FitMi is a full-body exercise tool that transforms classic therapy into a fun, interactive experience.

How Does it Work?

With FitMi, patients use two wireless “Pucks” to interact with therapeutic exercise apps on a Flint Tablet or a PC or Mac.

Each Puck contains 10 different sensors that track your patients’ movements in order to provide real-time feedback as they exercise in the apps.

As your patients exercise, RehabStudio provides real-time visual, auditory, and haptic feedback after each rep is completed, keeping your patients engaged throughout their therapy session.



CUSTOM REGIMEN

	REPS	HIGH SCORE
<input type="radio"/> Reach to Target #2	13	46
<input type="radio"/> Rotation	21	50
<input type="radio"/> Shoulder Flexion	11	45
<input checked="" type="radio"/> Shoulder Rotation	30	30
<input checked="" type="radio"/> Twists	20	20
<input type="radio"/> Standard Crunch	19	19
<input type="radio"/> Wrist Flexion	30	36
<input type="radio"/> Key Pinch	61	61
<input type="radio"/> Finger Tapping	25	25
<input checked="" type="radio"/> Twisting	120	120

At the end of each session, FitMi displays the number of reps your patients completed correctly for each exercise, and it lets them track their performance over time.