



*FITMI WITH 10" TABLET  
(CLINICS ONLY)*



*FITMI WITH 21" MONITOR*

## PRODUCT DETAILS

### **What is FitMi?**

FitMi is a new rehab tool that makes tabletop and lower extremity exercise much more exciting for your patients.

The device comes with two small wireless pucks that connect to a proprietary tablet. Patients then use the two pucks to complete exercises shown in the app. The therapist can tailor the difficulty of each exercise to the patient's motor impairment.

### **What is included with FitMi?**

Each FitMi station comes with:

- 2 FitMi Pucks
- 1 FitMi charging dock
- 1 FitMi Receiver
- Choice of 10" Tablet or 21" Monitor\* (\*Clinic-use only)

### Who is FitMi intended for?

FitMi is designed for all levels of impairment after a neurologic injury such as:

- ✓ Stroke
- ✓ Spinal Cord Injury
- ✓ Cerebral Palsy
- ✓ Traumatic Brain Injury
- ✓ Neurologic & Muscular Injury
- ✓ Developmental Disability

### How does the treatment work?

As your patients exercise, FitMi provides real-time visual, auditory, and haptic feedback when each rep is completed successfully. This interactive feedback will keep your patients engaged throughout their regimen.

Then, at the end of each session, the app displays the number of reps your patients completed correctly for each exercise and lets them see how they stack up against their best score.

### How much movement is needed?

FitMi is designed for all impairment levels. Even individuals with complete paralysis can benefit from using FitMi passively (meaning, you assist your affected limb), which stimulates the brain.

### Is any special training required?

FitMi is ready to use right out of the box with no special training required. The entire unit takes less than 5 minutes to set up and is designed to be user and clinician friendly. If you need assistance, you can contact Flint Rehab Devices at (800) 593-5468

- Flint's instructional videos can be found here: <https://www.flintrehab.com/howto/>

### What are the benefits of the Flint Tablet & Flint Monitor?

Flint offers dedicated tablets & monitors which connect with all of their devices. Each workstation has the app software preinstalled on it so that it works right out of the box and complies with FDA regulations. You simply push the power button on the tablet to start the therapy game. It does not need an internet or WiFi connection to work.

**I want to order both FitMi and MusicGlove. Does that mean I have to buy two tablets?**

No. If you would like both FitMi and MusicGlove, you can 'Add On' one device without a tablet. All Flint tablets come with both FitMi and MusicGlove software preinstalled.

For example, if you just purchase a FitMi, it will come with both FitMi and MusicGlove software installed so that you're ready to go if you decide to get a MusicGlove in the future.

## CLINIC INFORMATION

**What's the difference between the home and clinic software for FitMi?**

In the home version of FitMi, exercises are unlocked when patients are ready for them using our Adaptive Challenge Algorithm.

In the clinic version of FitMi, clinicians can create custom regimens using a drag-and-drop builder; and they can create unlimited profiles for each patient.

**What is the recommended regimen for FitMi?**

- For clinicians, we recommend the exercises in the sequence they are displayed in the quick play mode. Clinicians have the ability of creating their own regimen.
- For patients, we recommend starting with the initial 3 exercises of each body part. You may alternate the body parts by day. We suggest a 30 min workout session on a daily basis. If the users can't tolerate a complete 30 min session, they can do 15 min in the morning and then 15 min in the afternoon.

**How is the FitMi cleaned?**

You can use a standard disinfectant cloth to clean the outside of the pucks. To clean the charging dock, make sure it is disconnected from an outlet. Then use a standard disinfectant cloth to wipe the docking station.

**CPT Codes (Current Procedural Terminology)**

97112 - Neuromuscular reeducation

97110 - Therapeutic procedure/exercise

## PRODUCT SPECIFICATIONS

### **FDA Listed Classification**

Class 1; Exerciser; Non-Measuring (Regulation Number: 890.5370)

Direct Access: <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfrl/details.cfm?lid=419880>

### **UNSPSC (*United Nations Standard Products and Services Code*)**

1. 42251500 – Cognitive and dexterity and perceptual and sensory evaluation and therapy products (Specifically: 42251503 – Therapeutic games)
2. 42251600 – Rehabilitation exercise devices and equipment **HCPCS**

#### **CODE:**

- A9270 - Non-Covered Item or Service
- E1399 – “FitMi Therapy” (use for Veterans Health Administration)

### **What is FitMi made of?**

FitMi is made of compliant medical grade ABS durable plastic.

### **Where is FitMi manufactured?**

Taiwan

## MANUFACTURER INFORMATION

### **Flint Rehabilitation Devices, LLC**

18023 Sky Park Circle, Suite H2, Irvine CA 92614

Tel: 1(800) 593-5468 / Fax: (949)419-1394 [www.flintrehab.com](http://www.flintrehab.com)

#### SALES

Devon Gerwolls  
sales@flintrehab.com  
(949) 667-0140 EXT 1

#### SUPPORT

Daniel Monroy  
support@flintrehab.com  
(949) 667-0140 EXT 2

#### MARKETING/MEDIA

Kari Dahlgren  
info@flintrehab.com  
(949) 667-014



## WARRANTY INFORMATION

### **Abridged-version of Limited Warranty**

Flint warranties all components of each product when purchased new and unused to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. This warranty includes labor, materials, and shipping charges incurred in replacement or repair of the product.

### **Full-version of Limited Warranty**

This warranty is extended only to the original purchaser who purchases this product when new and unused from Flint Rehabilitation Devices, LLC (hereafter, Flint) or a dealer. This warranty is not extended to any other person or entity and is not transferable or assignable to any subsequent purchaser or owner. Coverage under this warranty will end upon any such subsequent sale or other transfer of title to any other person.

This warranty gives you specific legal rights and you may also have other legal rights which vary from state to state.

**Flint warrants all components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase from Flint or a dealer**, with a copy of the seller's invoice required for coverage under this warranty. If within such warranty periods any such product shall be proven to be defective, such product shall be repaired or replaced, at Flint's option.

This warranty does not include any shipping charges incurred in replacement part installation or repair of any such product. Flint's sole obligation and your exclusive remedy under this warranty shall be limited to such repair and/or replacement. For warranty service, please contact the dealer from whom you purchased your Flint product. In the event you do not receive satisfactory warranty service, please write directly to Flint. Provide dealer's name, address, the product model number, date of purchase, indicate nature of the defect and, if the product is serialized, indicate the serial number. Do not return products to our factory without our prior consent.

### **Extended Warranty**

The purchaser may purchase an additional one (1) year warranty for an additional 15% of the total product price.