

About Our Company

We make tools to spark recovery

At Flint Rehab, we believe rehabilitation is about more than just a regimen. It's about rebuilding a life full of the activities you love – and having fun while you're at it! That's why we've dedicated ourselves to using the latest clinical research and innovative ideas to create effective, engaging tools to spark recovery.



Flint Rehabilitation Devices LLC

18023 Sky Park Circle Ste. H2, Irvine CA 92614

www.FlintRehab.com

Office: 1(800) 593-5468 Fax: (949) 419-1394

Email: Support@FlintRehab.com

We Make It Easy on You

Our practical design makes it...

Easy to Request

Our products can be requested for direct veteran use or purchased for clinical treatment use as therapy equipment

Easy to Use

Everything come ready to go, right out of the box. No special training required!

Easy to Improve

Our tools make therapy fun and effective, which motivates you to achieve a higher recovery.

Total Continuum of Care

Versatile rehab tools for all levels of patient care

From inpatient to outpatient to home health/veteran use, our tools will provide an individual, clinician, and therapy department a seamless transition of care.



Inpatient /
Specialized Hospital



Outpatient /
Community-Based Clinic



Home Health / Veteran Use

MUSICGLOVE HAND THERAPY

Music-Based Hand Therapy Makes Rehab Fun

MusicGlove is an FDA listed device that includes a wearable sensorized glove and a therapy-based musical game on an included tablet.



MusicGlove motivates patients to perform hundreds of functional hand and finger repetitions while playing an engaging musical game. It is gradeable to the patient's level of impairment and tracks their progress after each therapy session. The device is safe, easy to understand, and can be used in the clinic with a therapist or at home without supervision.

How Does It Work?

Step 1: Select your grips

Choose which fingers you'd like to work on. Lower level patients can start with one grip.

Step 2: Select your difficulty level

Choose your difficulty level, then choose how long you'd like to workout for.

Step 3: Play!

Notes will float down the screen and you will attempt to make the correct pinching movement as each note hits the bottom of the screen.







FOR PATIENT USE



Includes:

- 1 MusicGlove (see the following page to find your size)
 - 10" touchscreen tablet with pre-installed software for 1 user
 - Headphones and connector cable
 - 1-year service contract (coverage on all software, hardware, and firmware)

*Eligible for coverage. See page 5 for details


MUSICGLOVE HAND THERAPY

Clinical Trials

Exercise with MusicGlove has been clinically proven to:

- Improve hand function in 2 weeks when used for 3 hours per week
- Lead to functional gains such as opening a door, washing dishes, typing, and using the restroom independently
- Motivate safe, high-intensity movements that initiate neuroplasticity in the brain

Subjects sustained these improvements at a one month follow-up and found the MusicGlove more motivating than other therapies.



FOR CLINIC USE

MusicGlove Clinic Portable Suite with 10" Tablet



Includes:

- 10" touchscreen tablet with pre-installed software
 - 6 MusicGloves (S, M, L pairs)
 - Software for unlimited users
 - Headphones and connector cable
- 1-year service contract (coverage on all software, hardware, and firmware)

FOR CLINIC USE

MusicGlove Clinic Stationary Suite with 21" Workstation



Includes:

- 21" touchscreen workstation with pre-installed software
 - 6 MusicGloves (S, M, L pairs)
 - Software for unlimited users
- Headphones and connector cable
- 1-year service contract (coverage on all software, hardware, and firmware)

MUSICGLOVE SIZING CHART & ELIGIBILITY REQUIREMENTS


Is MusicGlove Covered for Patients?

Yes! The home version of MusicGlove is eligible for coverage through the Prosthetic and Sensory Aids Service if a VA therapist believes there is proper medical justification for the device.

Eligibility Requirements

To use MusicGlove hand therapy, you need the ability to touch your thumb to at least one of your fingertips or side of your index finger.






Sizing Chart Instructions

If you have a printer available, print this chart out and **place your index finger** over each sizing model to determine your best fit. Make sure that your finger does NOT go outside the black lines. If it does, then go up a size.

If a printer is not available, then compare the tip of your **index finger** to the following standard American coins.



FITMI SMART NEURO REHAB DEVICE

FitMi Full-Body Therapy Is Like Wii for Rehab

FitMi is a rehabilitation tool that contains two sensorized pucks and a therapy based app that makes tabletop exercise more engaging for patients. Therapists can schedule 40 different exercises for patients like reach-to-target, rolling, stacking, and gripping.

The FitMi app logs all patient data, and provides real time visual and haptic feedback for the patient after each repetition is complete. The device is safe, easy to understand, and can be used in the clinic with a therapist or at home without supervision.

How Does It Work?

BEGIN FlintRehab

Step 1: Create Your Regimen

Therapist select exercises for

patients for hands, arms, core

or legs.



Step 2: Start Exercising

Using the blue and yellow pucks, complete the exercise shown on the screen.

Step 3: Continue to Next Exercise

After completing the # of reps prescribed, patients progress to the next exercise selected by the therapist.

FOR PATIENT USE

FitMi Home Portable Suite with 10" Tablet



Includes:

- 2 FitMi Pucks
- 10" touchscreen tablet with preinstalled software for patient use
 - 1 FitMi charging dock and Reciever
- 1-year service contract (software upgrades, replacement of any damaged


FITMI SMART NEURO REHAB DEVICE

Results through Therapeutic Repetition

Studies show that the amount of repetition completed during normal therapy is very low — in the 30's-40's.

FitMi makes it fun and easy to complete a high number of repetitions — which helps activate neuroplasticity.

In our initial pilot studies, we found that even individuals with severely limited movement were able to complete over 400 exercise repetitions in a 30-minute session.



No Eligibility Requirements

FitMi is adjustable to all stages of recovery following neurological injury. You can also adjust the difficulty of each exercise so that your patients are always optimally challenged.

FOR CLINIC USE FitMi Clinic Portable Suite with 10" Tablet Includes: 2 FitMi Pucks • 10" touchscreen tablet with preinstalled software for clinic use • 1 FitMi charging dock and Reciever • 1 USB drive with RehabStudio clinic software • 1-year service contract (software upgrades, replacement of any damaged hardware)

FOR CLINIC USE FitMi Clinic Stationary Suite with 21" Workstation CUSTOMIZE YOUR WORKOUT Includes: 2 FitMi Pucks 21" touchscreen workstation with preinstalled software for clinic use • 1 FitMi charging dock and Reciever • 1 USB drive with RehabStudio clinic software • 1-year service contract (software upgrades, replacement of any damaged hardware)

MUSICAL VARIATION AND ACCESSORIES



Song Extension Pack for MusicGlove - Pop

The Pop Song Pack adds 12 new pop songs to your MusicGlove experience. This adds twice as much variety to your regimen!

MusicTouch App for FitMi

MusicTouch is a music-based therapy app for FitMi. With this app, you get to play along to a musical game while you use your FitMi pucks to 'hit' each note that floats down the screen. Our users have a lot of fun playing MusicTouch!

*Tablet not included





Flint Travel Bag

Want a convenient way to transport your rehabilitation tools? Our Flint bag is perfect for carrying your rehabilitation tools from home to the clinic. It also makes a great travel tool, so that you can take your home therapy with you wherever you go.

FAQ & WARRANTY

Frequently Asked Questions

What's the difference between MusicGlove and FitMi?

MusicGlove is best for hand rehabilitation and FitMi is best for full-body rehabilitation.

Can I use MusicGlove with no hand movement?

MusicGlove is intended for patients who can touch their thumb to at least one of their fingertips. However, if you cannot make that movement but are still interested in MusicGlove, then you can try using the device passively.

Is any special training required?

FitMi and MusicGlove both come ready for use right out of the box with no special training required. They both take less than 5 minutes to set up and are designed to be user friendly. If you would like help setting up, simply call us at 1(800) 593-5468.

What's the difference between the home and clinic software for FitMi?

In the home version of FitMi, exercises are unlocked when patients are ready for them using our Adaptive Challenge Algorithm. In the clinic version of FitMi, clinicians can create custom regimens using a drag-and-drop builder; and they can create unlimited profiles for each patient.

I want to order both FitMi and MusicGlove. Does that mean I have to buy two tablets?

No. If you would like both FitMi and MusicGlove, you can Add On one device without a tablet. All Flint tablets come with both FitMi and MusicGlove software preinstalled. For example, if you just purchase a FitMi, it will come with both FitMi and MusicGlove software installed so that you're ready to go if you decide to get a MusicGlove in the future.

Limited Warranty

This warranty is extended only to the original purchaser who purchases this product when new and unused from Flint or a dealer. The warranty is not extended to any other person or entity and is not transferable or assignable to any subsequent purchaser or owner.

Coverage under this warranty will end upon any such subsequent sale or other transfer of title to any other person. This warranty give you specific legal rights and you may also have other legal rights which vary from state to state.

Flint warranties all components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase from Flint or a dealer, with a copy of the seller's invoice required for coverage under this warranty.

The comprehensive warranty covers all components of the system including hardware, so ware, and firmware. If within such warranty periods any such product shall be proven to be defective, such product shall be repaired or replaced, at Flint's op on. This warranty includes labor, materials, and shipping charges incurred in replacement or repair of the product. Flint's sole obligation and your exclusive remedy under this warranty shall be limited to such repair and/or replacement.

HOW TO ORDER

Thank you for your interest in our rehabilitation tools!

To purchase our products, please contact one of our premier government partners:



Service Disabled Veteran Owned
Small Business

Pisces Healthcare Solutions

Address: 2670 SW 76th Ave Miami, FL 33155

Contract Vehicles: FSS #: V797D-30191

Order Information: customerservice@pisceshealth.com Phone: (888) 505-3627 Fax: (888) 505-1987



Woman Owned Small Business

Tryco Inc.

Address: 6736 Old McLean Village Dr. McLean, VA 22101

Contract Vehicles:
FSS #: V797D-70077

DAPA Equipment: SP0200-09-H-0050

DAPA Med/Surg: SP0200-02-H-0042

ECAT Equipment: SPE2DH-16-D-0017

Order Information: tryco@tryco.org Phone: 800-934-3452 Fax: 703-734-0511

Manufacturer Information

Flint Rehabilitation Devices, LLC 18023 Sky Park Circle Suite H2

18023 Sky Park Circle Suite H2 Irvine, CA 92614

1-800-593-5468 support@flintrehab.com

Thank you for your interest in our rehabilitation tools.

For further information please visit us at: FlintRehab.com

You can also contact a Flint representative

Office: 1(800) 593-5468

Fax: (949) 419-1394

email: support@FlintRehab.com





Flint Rehabilitation Devices LLC 18023 Sky Park Circle Ste. H2 , Irvine CA 92614